

TD5B Oil Test Centrifuge



Features:

1. Programmable microprocessor control system, big LCD indicates rotor No., speed, time, temperature, RCF, rate of speed as well as fault information.
2. Adopts band heater, fast heating speed, temperature uniform, with the function of temperature-control and constant.
3. Centrifugal tube made in special skill, anti-broken when work.
4. Great torque brushless motor, maintenance-free, no powder, fast acceleration and deceleration, low noise, which improve the stability a lot.
5. Electric lid lock, over-speed and over-temperature protection
6. 3 phase to keep the best performance.
7. Stainless steel chamber, safe and reliable

Technical Parameters:

Max speed	4000rpm	Max RCF	2498xg
Max volume	4x100ml	Noise	≤60dBA
Timer	1~99min	Net weight	67KG
Dimension(LxWxH)mm	685×500×405mm	Temperature Range	Room temperature:+10~80℃
Speed accuracy	±20rpm	Package	Wooden box

