

TGM12 Hematocrit Centrifuge



Application:

TGM12 is specific centrifuge for Dengue Fever..


Features:

1. Adopts micro brushless motor, free-maintenance, no powder pollution, micro-computer control, LCD displays the time and speed.
2. Electric lid lock, easy and reliable.
3. Over-speed and imbalance protection.
4. Frame is 3 tires protection steel jacket, with stainless steel chamber.

TGM12 Technical Parameter:

Max. Speed	12000rpm
Max. RCF	15800
Max. Capacity	24pcs Capillary
Time Range	0~99min59s
RPM/RCF Convert	Yes
Noise (dB)	≤ 55
Temperature	Normal
Acc/Dec	10 Kinds
Speed Accuracy	± 20r/min
Temperature Accuracy	/
Voltage(V/Hz)	AC 220V/110V 50HZ/60HZ
Size (W x D x Hmm)	355×270×205mm
Net Weight(Kg)	13KG
Certificates	CE,ISO & Calibration report are available

Matched Rotors for TGM12

Rotor Photo	Max speed(rpm)	Max capacity	Max RCF (×g)
	12000	24 pieces vessel	15800