

## FYQ-100 Gel Card Incubator

### Application

Digital incubator for the Gel cards. Blood card/ID card incubator holds 2 independent incubation zones with a capacity for 12 Gel cards, It can be suitable for many types of gel cards from different manufacturers.



### Features

1. Keep 24-Gel Cards in temperature of 37 degree.
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital display to set up time and keep the temperature
4. Audible programmable alarm to notify end of incubation period

### Technical Parameter

Ambient Temperature	10~37℃
Temperature uniformity	≤±0.5℃
Max. Capacity	24 Blood cards
Size (W x D x Hmm)	390×280×140mm
Voltage(V/Hz)	AC 220V/110V 50HZ/60HZ
Incubation Temperature	37℃±0.1℃
Temperature Accuracy	<0.5℃
Rated Power	250W
Net Weight	6.2KG
Time Range	0~99min59sec
Certificates	CE,ISO & Calibration report are available